

Analytical Methods for Detection and Quantitation of Explosives Residues

Robert P. Jones

Chemistry Team Leader

Environmental Chemistry Branch

Environmental Laboratory, ERDC, Vicksburg, MS




US Army Corps
of Engineers

Engineer Research & Development
Center

Analytical Techniques

- High Performance Liquid Chromatography with UV detection (HPLC-UVD)
- Gas Chromatography with Electron Capture Detection (GC-ECD)
- Liquid Chromatography coupled with Mass Spectrometry (LC/MS)
- Gas Chromatography with Mass Selective Detection (GC/MSD)




Extraction Techniques

- Sonication - soils and sediments
- Direct injection - aqueous samples
- Salt-out extraction - aqueous samples
- Solid Phase Extraction - aqueous samples
- Solid Phase Microextraction (SPME) – aqueous samples




Analytical Methods

- SW-846 Method 8330 (current “gold standard”) uses dual column reverse phase HPLC with UV detection.
- SW-846 Method 8095 (Update IV) uses dual column high resolution GC with electron capture detection.
- Development work underway with LC/MS using Method 8321
- Possible role for GC/MS using modified Method 8270



SW-846 Methods

- Method 8330 - Nitroaromatics & nitramines by HPLC
- Method 8331 – Tetrazene by HPLC
- Method 8332 – Nitroglycerine by HPLC
- Method 8095 – Explosives by GC
- Method 8091 – Nitroaromatic & cyclic ketones by GC
- Method 8321 – Solvent extractable nonvolatile compounds by HPLC/TS/MS or HPLC/UVD
- Method 3535 – Solid Phase Extraction



Method 8330

- Current “gold standard” developed by Tom Jenkins *et al.* at CRREL in late 1980’s
- Dual column reverse phase (C-18 & CN) HPLC with UV detection
- Broad range of nitroaromatics and nitramines
- Direct injection of aqueous samples
- HPLC accurately quantitates extract concentrations as low as 20 ppb
- Aqueous DLs ~20 ppt when coupled with SPE (Method 3535)




Explosives Currently Analyzed by HPLC at ECB

- Method 8330 targets: HMX; RDX; 1,3,5-TNB; 1,3-DNB; Tetryl; NB; TNT; 4ADNT; 2ADNT; 2,4-DNT; 2,6-DNT; 2-NT; 3-NT; 4-NT
- PETN; Picric Acid; Nitroglycerine; CL-20
- TNT metabolites Azoxytoluenes, Diaminonitrotoluenes & Triaminotoluenes
- Nitroso metabolites of RDX (MNX,DNX,TNX) & HMX (MN-HMX)




Method 8095

- Relatively new developmental standard method appears in Update IV of SW-846
- Dual capillary column GC with electron capture detection
- Instrument accurately quantitates extract concentrations as low 1 ppb for some targets
- Aqueous DLs ~1 ppt when coupled with SPE (Method 3535)
- Incorporates a surrogate; ECB uses 3,4 DNT
- Targets: HMX; RDX; 1,3,5-TNB; 1,3-DNB; Tetryl; NB; TNT; 4ADNT; 2ADNT; 2,4-DNT; 2,6-DNT; 2-NT; 3-NT; 4-NT; 3,5-DNA; PETN; NG





Method 8095 Technical Issues



- Thermally labile explosives don't like hot GC injection ports and ovens
- HMX is very unstable; currently only one analytical column gives reproducible chromatography
- PETN also problematic giving poor chromatography with possible thermal degradation
- Cleanup techniques for extracts need to be developed
- Need frequent cleaning of GC injection port to minimize degradation
- Stability of low-level explosives concentrations in sample extracts



Environmental Chemistry Branch
Explosives by GC/ECD
(Restek Rtx-TNT Capillary Column)




Environmental Chemistry Branch
Explosives by GC/ECD
(Restek Rtx-TNT2 Capillary Column)





Explosives Analysis By GC/MSD

- Potential advantages
 - identification of unknowns when in SCAN mode
 - single analytical column
 - low detection limits (~ 100 ppb) in SIM mode
- Disadvantages
 - degradation of thermally labile compounds (HMX)
 - high flow rates present problems with MS vacuum
 - relatively high detection limits in SCAN mode



Environmental Chemistry Branch
Explosives by GC/MS
Total Ion Chromatogram - HP-5




US Army Corps
of Engineers

Engineer Research & Development
Center

Explosives by LC/MS

- For routine analytical work Method 8321 can be used for method development
- Mixtures separated into components by LC providing retention time data prior to introduction into MS
- Both SCAN and SIM modes can be used
- Can be coupled with UV detection
- Low detection limits
- Identification of unknowns compounds for which analytical standards are not commercially available




Explosives by LC/MS

- ECB's ThermoFinnigan system includes
 - LC/MS - LCQ Deca
 - Autosampler - Spectrasystem AS3000
 - LC Detector – Spectrasystem UV6000
 - Software – Xcalibur data system
- MS detector can operate in ESI (Method 8321) or APCI mode
- MS/MS analysis is two-stage process identifying characteristic fragment ion upon activation of isolated precursor ion
 - allows identification of targets in complex matrices
- Sample can be introduced by direct infusion without using LC for separation; mainly for research applications




Explosives by LC/MS

- Probably not a viable solution for routine analytical
- Relatively expensive
- Can be useful with problem samples
- Good tool for research
- Valuable for qualitative analysis of unknowns
- Can be used for confirmation when not possible by HPLC alone




Detection Limits

- Regulatory limits or project specific requirements dictate detection limits needed
- Influenced by sample matrix, nature of target compound, sample prep & analytical techniques
- Sample prep: Concentration techniques
→ 1L water sample → 10 mL SPE extract → 20 uL microconcentrate = 50,000:1
- Analytical technique sensitivity:
GC/ECD > HPLC/UVD > GC/MS or LC/MS (SIM) > GC/MS or LC/MS (SCAN)



Summary


- **HPLC** → mature methodology; good sensitivity; ECB's broadest range of targets
- **GC** → developing technology; currently offers highest sensitivity; techniques need to be honed
- **GC/MSD** → offers some advantages to GC & HPLC; probably not viable for routine analysis of explosives
- **LC/MS** → useful for problem samples & research applications; sensitive & selective; probably not viable for routine analysis of explosives



Environmental Chemistry Branch

Mission

- Quality assurance: split sample quality assurance analysis, technical assistance, and problem solving to support USACE environmental HTRW programs
- Analytical chemistry support: water quality analysis and in-house HTRW project analysis
- Research: DOD relevant topics
- Research support: chemical analysis and chemistry technical support for ERDC research




Environmental Chemistry Branch Resources

- Personnel: 31 Federal Government employees, 4 contract employees, 14 contract students
- Branch Chief: Doug Taggart
- Four team leaders: Dave Splichal, Prem Arora, Laura Percifield, and Bobby Jones
- Equipment: all major routine analytical equipment plus specialized equipment such as ICP/MS, LC/MS
- Excellent analytical chemistry facilities in two locations




Standard Analytical Chemistry Capabilities

- Volatile and semivolatile organics by GC/MS
- Pesticides and herbicides by GC
- Metals by ICP, ICP/MS, and GFAA
- Nutrients and anions by ion chromatography and autoanalyzers
- Explosives by HPLC
- Petroleum hydrocarbons by GC
- Polychlorinated biphenyls (PCBs) by GC



Specialized Analytical Capabilities

- Metals by ICP/MS
- PCB congeners
- Low-level PAHs by GC/MS (SIM)
- Explosives (TNT and RDX) breakdown products analysis by HPLC
- White phosphorus
- Metals speciation (As, Fe, Cr, Hg)
- Other methods under development: LC/MS of explosives, perchlorate by IC, methyl mercury, explosives by GC/ECD, tributyl tin by GC



Contacts

- Dr. Douglas B. Taggart, Branch Chief, (402) 444-4300; e-mail: Douglas.B.Taggart@usace.army.mil
- Robert P. Jones, team leader, (601) 634-4098 or 1-800-522-6937 (voice); (601)-634-2742 (fax); e-mail: Robert.P.Jones@erdc.usace.army.mil
- Laura J. Percifield, team leader, (402) 444-4314; e-mail: Laura.J.Percifield@nwo02.usace.army.mil
- Prem N. Arora, team leader, (402) 444-4318; e-mail: Prem.N.Arora@nwo02.usace.army.mil
- Dave Splichal, team leader, (402) 444-4318; email: David.E.Splichal@nwo02.usace.army.mil

